Lugol's Solution with 5 % Iodine (Concentrated)

Principle

Lugol's solution is an aqueous solution of iodine and potassium iodide. It acts as a mordant in the Gram-Hucker staining method. A mordant is a substance that enables or enhances the staining of a biological sample by a dye. In Gram-Hucker staining, the dye is usually crystal violet, and the role of iodine as a mordant is essential for the proper performance of the staining procedure.

At a higher concentration (5 %), it is widely used in protozoology, from rapid fixation to general staining of human intestinal protozoa. It also reveals starch, glycogen, and other polysaccharides.

This concentration is used as a stock solution to prepare, by dilution, the Lugol's solution used in Gram staining. This staining method is qualitative.

Material

- Biological samples (sputum, pus, cultured bacterial cells).

Reagents

Product code	Product name
252532	Gram-Hucker's Crystal Violet Oxalate solution (CE-IVD) for clinical diagnostics (*)
251803	Alcohol-Acetone 7:3 (CE-IVD) for clinical diagnostics (*)
252531	Gram-Hucker's Safranin O solution (CE-IVD) for clinical diagnostics (*)

Procedure

In Gram staining, biological samples (sputum, pus, cultured bacterial cells) are used, usually prepared as smears and fixed with methanol or by flame fixation. After staining with crystal violet, the preparation is treated with a 0.4 % Lugol's solution so that it penetrates the cell wall and interacts with the crystal violet previously applied. The preparation is then decolorized and safranin is used as a counterstain. The preparation can then be observed under a light microscope, revealing Gram-positive cells (violet-blue) and Gram-negative cells (red-pink).

The procedure for Gram staining is as follows:

- Heat-fix the microbial smear, taking care not to burn the preparation. Apply short passes while monitoring the slide temperature by touching the back of the hand.
- Cover the preparation with Crystal Violet Oxalate Solution according to Gram-Hucker (252532) for 1 minute.
- Rinse with running water.
- Cover the preparation with the 0.4 % diluted working solution for 1 minute.
- Rinse again with water.
- Decolorize drop by drop for no longer than 1 minute with Alcohol-Acetone 7:3 (251803).
- Rinse again with water.
- Cover the preparation with the safranin solution.
- Rinse with water and allow to dry before microscopic examination.

Results

Gram-positive bacteria appear deep violet. Gram-negative bacteria appear pink. False results may occur due to handling errors or unusual characteristics of specific organisms. Particular care should be taken during fixation and decolorization, since excessive treatment may damage the cell wall and lead to excessive removal of crystal violet, causing false-negative results.

Technical note

The microscope used should meet the requirements of a clinical diagnostic laboratory. If an automatic staining instrument is used, the manufacturer's operating instructions and software instructions must be followed.

Sample preparation

All samples must be handled according to the current state of technology. All samples must be clearly and unambiguously labeled.

Diagnosis

Diagnoses must be established only by authorized and qualified personnel. Appropriate controls should be included in every application to rule out erroneous results.

Storage

Keep in tightly closed containers protected from light. Store at room temperature (+15 °C to +25 °C).

Expiry

The product, when stored at the indicated temperature and in a tightly closed container, may be used until the expiry date stated on the packaging.

Notes on use

To avoid errors, staining must be carried out by trained personnel. For professional use only. National regulations on occupational safety and quality assurance must be observed. Consult the Safety Data Sheet before use.

Waste disposal information

Used and expired solutions must be disposed of as hazardous waste in accordance with local regulations. Questions regarding disposal or procedures may be sent by e-mail to: info.es@itwreagents.com.

Within the EU, the applicable provisions are based on Regulation (EC) No. 1272/2008 of the European Parliament and of the Council of 16 December 2008 on the classification, labeling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC and amending Regulation (EC) No. 1907/2006.

Classification of hazardous substances

Refer to the hazardous substance classification on the label and the information provided in the Safety Data Sheet.

Revision history

Version March 2022	2022
Version August 2026	2026 Transferred to an updated template

Manufacturer

Panreac Química S.L.U. - An ITW Company
C / Garraf, 2 - Polígono Pla de la Bruguera
E-08211 Castellar del Vallès (Barcelona) Spain
Tel. (+34) 937 489 400 Fax (+34) 937 489 401

(*) Sanitary product for In Vitro Diagnostics in accordance with Regulation (EU) 2017/746 (IVDR)

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