

## RNase-ExitusPlus™

*RNase decontamination solution*

**Product codes A7153**

---

RNase-ExitusPlus™ is a fast and effective ready-to-use agent for laboratory surfaces and equipment. The decontamination starts instantly after spraying on a contaminated surface. The benefit of the product is, that it is neither toxic for humans nor corrosive for laboratory material. Competitor products use unhealthy and corrosive agents. RNase-ExitusPlus™ is a non-alkaline and non-carcinogenic cleansing solution that is highly active against RNase contamination. RNase-ExitusPlus™ has been demonstrated to inactivate more than 20 µg of RNase A dried onto the bottom of a microcentrifuge tube. RNaseExitusPlus™ is stable for approximately 12 months and heat resistant.

These are the new and unique characteristics of **RNase-ExitusPlus™**:

- 1) Catalytic and cooperative effects of the components cause a very rapid inactivation of protein and RNase molecules.
- 2) All components of RNase-ExitusPlus™ are readily biodegradable and not harmful or toxic to humans.
- 3) No aggressive mineral acids or alkaline substances are used. Equipment and materials are not damaged or corroded even after prolonged incubation times.
- 4) No toxic fumes. The reagent contains a low volume of alcohol only.
- 5) Elevated temperatures above approx. 50°C speed up the reaction and the efficiency / activity!

RNase-ExitusPlus™ is ready-to-use for eliminating RNase from any surface including the interior of microcentrifuge tubes. By following the simple decontamination instructions below, RNase is completely inactivated and removed. RNase-ExitusPlus™ should be stored at room temperature; at colder temperatures, a precipitate may form which is easily brought into solution at 37°C.

### Detailed instructions

**To decontaminate laboratory surfaces:** Apply RNase-ExitusPlus™ directly to the lab surface. Wipe thoroughly with a paper towel, rinse with water, and dry with a clean paper towel.

**To decontaminate laboratory apparatus:** Generously apply RNase-ExitusPlus™ to a paper towel and wipe all exposed surfaces of the apparatus thoroughly. Rinse with water and dry with a clean paper towel. To clean small parts, briefly soak them in RNase-ExitusPlus™, rinse with water and dry.

**To decontaminate plastic and glass vessels:** Add ample RNase-ExitusPlus™ to enable coating the entire surface of the vessel by swirling or vortexing. Discard the solution and rinse the vessels thoroughly two times with distilled water.

**To decontaminate pipettors:** Following the manufacturer's instructions; remove the shaft from the pipettor and remove seals and gaskets from the shaft. Soak the shaft for one minute in RNase-ExitusPlus™, rinse the shaft thoroughly with water, let dry and reassemble.

### Quality control

Aliquots of RNase A (10 µg) were dried down in reaction tubes for samples 1, 3, and 4 (**Fig. 1**). Afterward, RNase A samples were treated with 1 ml RNase-ExitusPlus™ (1) or H<sub>2</sub>O (3, 4) for 5 minutes at RT. Two washing steps with 1 ml of sterile water followed. Then 5 µg total RNA from *E. coli* was added into each tube. Into tube 4, a fresh aliquot of 10 µg of RNase A was added. All tubes were incubated for 30 min. at 37°C. Finally, a loading buffer was added, and samples were loaded onto a 1% agarose gel. As a control, 5 µg untreated total *E. coli* RNA (C) was included.




Fig. 1




Fig. 2

**Fig. 2. Analysis of autoclaved proteins without (-) and with (+) the addition of RNase-ExitusPlus™.** Test solutions of 10 mM Tris, pH 8.0 with BSA (bovine serum albumin) and RNase A were autoclaved at 120°C and 1.2 bar for 20 minutes after the addition of equal volumes of either sterile water (-) or RNase-ExitusPlus™ (+). Subsequently, aliquots of 10 µl with 1µg BSA or RNase A, respectively, were analyzed on a 4-12 % polyacrylamide gel and stained with Coomassie Brilliant Blue. The sample containing sterile water (-) shows no significant degradation of the proteins, while the addition of RNase-ExitusPlus™ (+) leads to almost complete degradation.