



PRODUCT CODE: 251563

Fehling's A Reagent for clinical diagnostics

for determination of glucose in urine and reducing sugars

TARIC 3822 90 00 00

SYNONYMS: Alkaline cupric solution

PHYSICAL DATA: Miscible with water D 1,024 •

INSTRUCTIONS FOR USE: 10 ml of solution A + 10 ml of solution B + 5 ml of solution of Potassium Ferrocyanide 10% + 25 ml of water, by heating, are reduced by 0.05 g of glucose

BIBLIOGRAPHY:

HAZARDOUS:



H: H318 • H410 •

P: P273 • P280 • P391 • P305+P351+P338 • P310 • P501 •

TRANSPORT REGULATIONS: UN: 3082 • ADR: 9/III • IMDG: 9/III • IATA: 9/III • PAX: 964 • CAO: 964 • (-) •

WEIGHT/VOLUME INFORMATION: 1l-1,024 kg 1kg-0,977 l

COMPOSITION:

Copper(II) Sulphate 5-hydrate	48,30 g
Sulphuric Acid 96%	1 ml
Water (s.q.m)	1 l

SPECIFICATIONS:

Identity :

Identity

passes test

Maximum limit of impurities

Suitability for determination of glucose

passes test