



**PRODUCT CODE: 253454**

**Weigert's Hematoxylin solution B for clinical diagnosis**  
for microscopy

---

TARIC 3203 00 10 00

**PHYSICAL DATA:** liquid, Clear, Greenish-yellow, Miscible with water • D 1,013 •

**INSTRUCTIONS FOR USE:** Mix equal parts of solution A and solution B. Stain with this Weigert's working solution for 10 minutes.

**BIBLIOGRAPHY:** Clark G., "Staining Procedures", 4th Ed., Williams and Wilkins ,Baltimore (1981) ; H. J. Conn , "Biological Stains", 9th Ed., Williamsand Wilkins , Baltimore , 609 (1977) ; R. García del Moral,"Laboratorio de Anatomía Patológica", 1ª Ed., McGraw-Hill , 133 (1993); S. Frankel , S. Reitman , A. C. Sonnenwirth , "Gradwohl#s clinicallaboratory methods and diagnosis" , 7th Ed., The C. V. Mosby Company 1, 943 (1970) •

**HAZARDOUS:**



H: H290 • H319 •

P: P234 • P264 • P280 • P305+P351+P338 • P406 • P501 • P337+P313 •

**WEIGHT/VOLUME INFORMATION:** 1l~1,013 kg      1kg~0,987 l

**COMPOSITION:**

|                                          |        |
|------------------------------------------|--------|
| Iron(III) Chloride 30% aqueous solution. | 4 ml   |
| Hydrochloric Acid 35%                    | 1 ml   |
| Water (s.q.m)                            | 100 ml |

**SPECIFICATIONS:**

Identity :  
Identity

passes test

Ed.: 3 . Vig.: 11.12.2015 .

Prod.: 253454