



PRODUCT CODE: 256879

Eosin Yellowish alcoholic solution 1% (CE-IVD) for clinical diagnostics

Solution for Hematoxylin-Eosin staining, in human and gynaecological samples., In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Red, Miscible with water • D 20/4 0,827 • Flash P.:20,8 °C(closed cup) •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 7.060 mg/kg •



H: H225 • H319 •

P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P280 • P303+P361+P353 • P370+P378 • P403+P235 • P264 • P337+P313 • P305+P351+P338 •

TRANSPORT REGULATIONS: UN: 1170 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-0,827 kg 1kg-1,209 l

COMPOSITION:

Eosin Yellowish	10 g
Acetic Acid glacial	1 ml
Water	50 ml
Ethanol 96 %	1000 ml

SPECIFICATIONS:

Identity :

Identity

passes test

Maximum limit of impurities

Suitability for hematoxylin - eosin stain

passes test