

PRODUCT CODE: 121347

Hexane, alkanes mixture for analysis

C₆H₁₄

M.= 86,18

CAS [92112-69-1]

EINECS 295-570-2

TARIC 2901 10 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Inmiscible with water. • D 20/4 0,67 • M.P.: -95 °C • B.P.: 60 - 70 °C • Flash P.: -22 °C • Ign. T.:240 °C • Vap. press. (20 °C) 177 hPa • Expl. limit0,6 %(V)7,7 %(V) •

BIBLIOGRAPHY: ACS XI •

HAZARDOUS: RTECS: MN 9275000 • LD50 oral rat 28.710 mg/kg • LC L0 mus 120000mg/m3 • VLA-ED 20 ppm72 mg/m3



H: H225 • H315 • H411 • H304 • H336 • H373 • H361f •

P: P210 • P233 • P240 • P241 • P242 • P243 • P261 • P264 • P271 • P273 • P280 • P301+P310 • P302+P352 • P303+P361+P353 • P304+P340 • P312 • P321 • P331 • P332+P313 • P362 • P370+P378 • P391 • P403+P233 • P403+P235 • P405 • P501 •

TRANSPORT REGULATIONS: UN: 1208 • ADR: 3(CONT.MARINO)/II • IMDG: 3(CONT.MARINO)/II • IATA: 3(CONT.MARINO)/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,67 kg 1kg~1,49 l

SPECIFICATIONS:

Minimum assay (as isomers C₆H₁₄) (G.C.) 95,0%

Maximum limit of impurities

Acidity 0,0003 meq/g

Non-volatile matter	0,001 %
Darkened substances by H2SO4	passes test
Sulfur compounds (as S)	0,005 %
Thiophene	passes test
Water (H2O)	0,01 %

Metals by ICP [in mg/Kg (ppm)]

Ag	0,05
Al	0,5
As	0,05
Au	0,05
B	0,02
Ba	0,1
Be	0,02
Bi	0,05
Ca	0,5
Cd	0,05
Co	0,02
Cr	0,02
Cu	0,02
Fe	0,1
Ga	0,02
Ge	0,05
Hg	0,05
In	0,05
K	0,1
Li	0,05
Mg	0,1
Mn	0,02
Mo	0,02
Na	0,5
Ni	0,02
P	5,0
Pb	0,1
Pt	0,02
S	0,2
Sb	0,02
Si	0,2
Sn	0,1
Sr	0,2
Ti	0,02
Tl	0,02
V	0,02
Zn	0,1
Zr	0,02