



PRODUCT CODE: 124860

Toluene / 2-Propanol / Water 100:99:1 v/v/v (Mixture TAN) (ASTM D 664) for analysis

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Colourless, Miscible with ethanol • D 20/4 0,82 • Flash P.:9 °C •

BIBLIOGRAPHY: ASTM Standard D 664-11a "Standard Test Method for Acid Number of Petroleum Products by Potentiometric Titration" •

HAZARDOUS:



H: H225 • H302 • H319 • H315 • H361d • H336 • H373 • H304 •

P: P280 • P301+P310 • P321 • P303+P361+P353 • P305+P351+P338 • P501 • P330 • P362+P364 • P405 •

TRANSPORT REGULATIONS: UN: 1993 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-0,822 kg 1kg-1,216 l

COMPOSITION:

| | |
|------------|--------|
| Toluene | 500 ml |
| 2-Propanol | 495 ml |
| Water | 5 ml |

SPECIFICATIONS:

Density at 20/4 0,820 - 0,830

Maximum limit of impurities

Acidity value 0,005

