

PRODUCT CODE: 125396

Phenol-1,1,2,2-Tetrachloroethane 60:40 w/w for analysis

for determination of the viscosity of polymers, according to ASTM D 4603, ISO 1628

TARIC 3822 90 00 00

SYNONYMS: 1,1,2,2-Tetrachloroethane-Phenol

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with ethanol • D 20/4 1,234 •

BIBLIOGRAPHY: Plastics-Determination of the viscosity of polymers in dilute solution using capillary viscometers ISO 1628-5:1998(E) • Plastics-Determination of the viscosity of polymers in dilute solution using capillary viscometers ISO 1628-5:1998(E), Standard Test Method for Determining Inherent Viscosity of Poly(Ethylene Terephthalate) (PET) by Glass Capillary Viscometer ASTM D4603-96 •

HAZARDOUS: VLA-ED (1,1,2,2-tetrac.) 1 ppm7 mg/m3 VLA-ED (phenol) 5 ppm20 mg/m3



H: H330 • H310 • H302 • H314 • H411 • H373 • H341 •

P: P260 • P262 • P264 • P270 • P271 • P501 • P273 • P280 • P284 • P301+P310 • P301+P330+P331 • P302+P350 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P310 • P320 • P321 • P322 • P330 • P338 • P361 • P363 • P391 • P403+P233 • P405 •

TRANSPORT REGULATIONS: UN: 2810 • ADR: 6.1/II • IMDG: 6.1/II • IATA: 6.1/II • PAX: 654 • CAO: 662 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-1,234 kg 1kg-0,81 l

SPECIFICATIONS:

Assay (as Phenol)	58-62 %
Assay (as 1,1,2,2-Tetrachloroethane)	38-42 %
Water (H ₂ O)	0,05 %

