



PRODUCT CODE: 126352

Amido Black 10B solution for the determination of protein content in milk for analysis

according to NF V 04-216

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Deep blue, Miscible with water • D 20/4 1,007 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l-1,007 kg 1kg-0,993 l

OBSERVATIONS: Storage between 15 and 25°C • Storage away from direct light. •

SPECIFICATIONS:

pH at 20°C	2,30-2,50
------------	-----------

ABS at λ 620 nm 1% sol. in H ₂ O	0,695-0,735
--	-------------