






PRODUCT CODE: 131159

Potassium Antimony(III) Tartrate 3-hydrate for analysis, ACS



M.= 667,87

CAS [28300-74-5]

EINECS 234-293-3

TARIC 2918 13 00 00

SYNONYMS: Antimony(III) Potassium Tartrate, Potassium Antimony(III) Oxide Tartrate, Potassium Antimonyl Tartrate

PHYSICAL DATA: Small crystals, White, Soluble in water 55 g/l at 20 °C D 2,6 • pH(50g/l)4,0 •

BIBLIOGRAPHY: Merck Index **12**, 741 13, 706 Safety **2**, **2895 B** • Römp **8**, **2012** • ACS **XI** • USP **41** •

HAZARDOUS: C.E: 051-003-00-9 • RTECS: CC 8825000 • LD L0 oral man 2 mg/kg • LD50 oral rat 115 mg/kg • LD50 ipr mus 45 mg/kg



H: H332 • H302 • H411 •

P: P261 • P264 • P270 • P271 • P273 • P501 • P301+P312 • P304+P340 • P312 • P330 • P391 •

TRANSPORT REGULATIONS: UN: 1551 • ADR: 6.1/III • IMDG: 6.1/III • IATA: 6.1/III • PAX: 670 • CAO: 677 • (E) •

SPECIFICATIONS:

Assay (Iodom.) 99,0-103,0%

Maximum limit of impurities

Acidity or alkalinity 0,02 meq/g

Insoluble matter in H₂O 0,02 %

Loss on drying at 105°C 2,7%

Chloride (Cl) 0,005%

Sulfate (SO₄) 0,01%

As 0,015 %

Ed.: 3 . Vig.: 15.01.2011 .

Prod.: 131159