





PRODUCT CODE: 142970

Sodium Stannate 3-hydrate pure



M.= 266,71

CAS [12209-98-2]

EINECS 234-724-5

TARIC 2841 90 85 90

PHYSICAL DATA: crystals, crystalline powder, White, Soluble in water 613 g/l at 16 °C pH(10% sol.)13 •

BIBLIOGRAPHY: Merck Index **12**, 8.826

HAZARDOUS: RTECS: JN 6345000 • LD50 oral rat 3.457 mg/kg •



H: H319 •

P: P264 • P280 • P305+P351+P338 • P337+P313 •

SPECIFICATIONS:

Assay (Sn)	42%
Alkalinity	29,0-32,5 %
Cu	0,001 %
Ni	0,001 %
Pb	0,001 %