



PRODUCT CODE: 145208

Pepsin 1:10.000 NF pure

CAS [9001-75-6] EINECS 232-629-3
TARIC 3507 90 90 90

PHYSICAL DATA: powder, White or slightly yellow, Soluble in water. Hygroscopic .

BIBLIOGRAPHY: Merck Index **13**, 7.225

HAZARDOUS:C.E: 647-008-00-6 .



H: H319 • H335 • H315 • H334 •

P: P261 • P264 • P271 • P280 • P285 • P302+P352 • P304+P340 • P304+P341 • P305+P351+P338 • P312 • P321 • P332+P313 • P337+P313 • P342+P311 • P362 • P403+P233 • P405 • P501 •

OBSERVATIONS: Storage between +2 and +8C, in a dry place, away from direct light. .

SPECIFICATIONS:

Proteolytic Activity	1:10000
Identity :	
Identity	IR passes test
pH of 5% solution	3 - 5

Maximum limit of impurities

Insoluble matter in water	passes test
Contamination Total aerobic microbial count (TAMC)	10000 cfu/g

Total combined yeast and moulds (TYMC)
Escherichia coli
Salmonella

100 cfu /g
absence/g
absence/10g

Ed.: 1 . Vig.: 10.05.2021 .

Prod.: 145208