



PRODUCT CODE: 15B496

Dimethylmethoxy Chromanol for synthesis

CAS [83923-51-7]

EINECS 893-588-0

TARIC 2909 30 38 90

PHYSICAL DATA: solid, crystalline, Off White to Brown, M.P.: 114 - 115 °C •

BIBLIOGRAPHY:

HAZARDOUS:



H: H317 • H411 •

P: P261 • P273 • P333+P313 • P280 • P501 • P321 •

TRANSPORT REGULATIONS: UN: 3077 • ADR: 9/III • IMDG: 9/III • IATA: 9/III • PAX: 956 • CAO: 956 • (-) •

SPECIFICATIONS:

Minimum assay	98%
Identity :	
Identity	IR passes test
Melting range	114-115°C
Residual solvents:	
2-Propanol	<5000 ppm
Toluene	<890 ppm
tert-Butyl Methyl Ether	<5000 ppm
Colour: Parameter L (D65)	70 - 95 *
Colour: Parameter C (D65)	8 - 18 *
Colour: Parameter h (D65)	75 - 100 *
Microbial Limits	
Total Aerobic Microbial Count (TAMC)	<100 cfu/g

Total combined yeast and moulds (TYMC)	<10 cfu/g
Heavy metals (as Pb)	<0,002 %

*** At the moment of the batch analysis**