

**PRODUCT CODE: 161347****Hexane, 95% alkanes mixture for synthesis** C_6H_{14}

M.= 86,18

CAS [92112-69-1]

EINECS 295-570-2

TARIC 2901 10 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Immiscible with water. • D 20/4 0,67 • M.P.: -95 °C • B.P.: 60 - 70 °C • Flash P.: -22 °C • Ign. T.: 240 °C • Vap. press. (20 °C) 177 hPa • Expl. limit 0,6 %(V) 7,7 %(V) •

BIBLIOGRAPHY: ACS XI •

HAZARDOUS: RTECS: MN 9275000 • LD50 oral rat 28.710 mg/kg • LC L0 mus 120000mg/m3 • VLA-ED 20 ppm/72 mg/m3



H: H225 • H315 • H411 • H304 • H336 • H373 • H361f •

P: P210 • P233 • P240 • P241 • P242 • P243 • P261 • P264 • P271 • P273 • P280 • P301+P310 • P302+P352 • P303+P361+P353 • P304+P340 • P312 • P321 • P331 • P332+P313 • P362 • P370+P378 • P391 • P403+P233 • P403+P235 • P405 • P501 •

TRANSPORT REGULATIONS: UN: 1208 • ADR: 3(CONT.MARINO)/II • IMDG: 3(CONT.MARINO)/II • IATA: 3(CONT.MARINO)/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-0,67 kg 1kg-1,49 l**SPECIFICATIONS:**

Minimum assay (as isomers C_6H_{14}) (G.C.)	95%
Non-volatile matter	0,001 %
Water (H ₂ O)	0,01 %

