





**PRODUCT CODE: 162325**

**2,4-Dinitrophenylhydrazine, 99% moistened with~ 33% of H2O for synthesis**



M.= 198,14

CAS [119-26-6]

EINECS 204-309-3

TARIC 2928 00 90 90

**SYNONYMS:** DNPH

**PHYSICAL DATA:** crystalline powder, Red, Insoluble in water Thermal decom. above 160 °C • M.P.: 200 °C •

**BIBLIOGRAPHY:** Merck Index **12**, 3.336 13, 3.312 Sax **DVC400** • Safety **2**, **1440 D** • Römp **8**, **973** • Beilstein **15**, **489 I**, **146 II**, **215 III**, **425 IV**, **380** • BRN 615586 • Fieser **1330 2176** •

**HAZARDOUS:** RTECS: MV 3325000 • LD50 ipr mus 450 mg/kg



H: EUH001 • H302 • H319 • H315 •

P: P264 • P270 • P280 • P301+P312 • P302+P352 • P501 • P305+P351+P338 • P321 • P330 • P332+P313 • P337+P313 • P362 •

**TRANSPORT REGULATIONS:** UN: 1325 • ADR: 4.1/III • IMDG: 4.1/III • IATA: 4.1/III • PAX: 446 • CAO: 448 • (E) •

**OBSERVATIONS:** Storage above 0°C •

**SPECIFICATIONS:**

Minimum assay (HPLC.)

99%

Identity :

Identity

IR passes test

Ed.: 3 . Vig.: 15.01.2011 .

Prod.: 162325