



PRODUCT CODE: 171569

Griess-Ilosvay's A Reagent

for detection of nitrites

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Colourless, Miscible with alcohol and water • D 20/4 1,045 •

INSTRUCTIONS FOR USE: Prepare the reagent by mixing equal volumes of A and B solutions. A pink colour develops in the presence of nitrites.

BIBLIOGRAPHY:

HAZARDOUS: CEE (CH₃COOH) 10 ml/m³25 mg/m³ VLA-EC (CH₃COOH) 15 ppm37 mg/m³ VLA-ED (CH₃COOH) 10 ppm25 mg/m³



H: H314 •

P: P260 • P264 • P280 • P301+P330+P331 • P303+P361+P353 • P501 • P304+P340 • P305+P351+P338 • P310 • P321 • P338 • P363 • P405 •

TRANSPORT REGULATIONS: UN: 1760 • ADR: 8/III • IMDG: 8/III • IATA: 8/III • PAX: 852 • CAO: 856 • (E) •

WEIGHT/VOLUME INFORMATION: 1l-1,045 kg 1kg-0,957 l

COMPOSITION:

Sulphanilic Acid	0,34 g
Acetic Acid glacial	30 ml
Water	64 ml

