



PRODUCT CODE: 173333

Vanadate-Molybdate Reagent
for determination of phosphates

TARIC 2807 00 00 00

PHYSICAL DATA: liquid, Clear, Yellowish, Miscible with water • D 20/4 1,149 •

BIBLIOGRAPHY: Standard Methods for the examination of water and waste water, 18th Ed., 4-113 (1992) •

HAZARDOUS: LD50 oral rat 2.140 mg/kg • TLV-TWA (H2SO4) 1 mg/m3 VLA-EC 3 mg/m3 VLA-ED (H2SO4) 1 mg/m3



H: H314 •

P: P260 • P264 • P280 • P301+P330+P331 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P310 • P321 • P338 • P363 • P405 • P501 •

TRANSPORT REGULATIONS: UN: 1760 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l~1,149 kg 1kg~0,87 l

OBSERVATIONS: Product controlled as a drug precursor. •

COMPOSITION:

Ammonium Molybdate 4-hydrate	2,97 g
Ammonium meta-Vanadate	0,15 g
Sulphuric Acid 96%	14,2 ml
Water (s.q.m)	100 ml

SPECIFICATIONS:

Suitability for determination of phosphates

passes test

Ed.: 2 . Vig.: 30.08.2017 .

Prod.: 173333