



**PRODUCT CODE: 175145**

### **OXI-OLEO-TEST**

for determination of the degree of deterioration of fried fats and oils, 60 tests

---

TARIC 3822 90 00 00

**PHYSICAL DATA:**

**BIBLIOGRAPHY:**

---

**HAZARDOUS:**



H: H225 • H332 • H302 • H314 • H318 • H336 •  
P: P210 • P233 • P240 • P241 • P242 • P243 • P260 • P261 • P264 • P270 • P271 • P280 • P301+P312 •  
P301+P330+P331 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P310 • P312 • P321 • P330 • P338 • P363 •  
P370+P378 • P403+P233 • P403+P235 • P405 • P501 •

---

**TRANSPORT REGULATIONS:** UN: 2924 • ADR: 3(8)/II • IMDG: 3(8)/II • IATA: 3(8)/II • PAX: 352 • CAO: 363 •  
(D/E) •

---

**SPECIFICATIONS:**

Comprised of:  
2x100 ml OXI-OLEO-TEST Reagent 1  
1x20 ml OXI-OLEO-TEST Reagent 2  
1 Plastic case  
2 Test tubes with stopper  
1 Scoop  
1 Colour scale  
1 Instructions sheet

