



PRODUCT CODE: 176131

Azidiol

for preserving milk samples

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Purple, Miscible with water • D 20/4 1,035 •

INSTRUCTIONS FOR USE: Dose 133 ul in 40 ml of milk or equivalent proportion.

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 27 mg/kg •



H: H301 • H413 •

P: P264 • P270 • P273 • P301+P310 • P321 • P501 • P330 • P405 •

TRANSPORT REGULATIONS: UN: 3287 • ADR: 6.1/III • IMDG: 6.1/III • IATA: 6.1/III • PAX: 655 • CAO: 663 • (E) •

WEIGHT/VOLUME INFORMATION: 1l~1,035 kg 1kg~0,966 l

COMPOSITION:

Chloramphenicol	0,75 g
Ethanol absolute	10 ml
Bromophenol Blue	0,35 g
Sodium Azide	18 g
tri-Sodium Citrate 5,5-hydrate	45 g
Water (s.q.m)	1000 ml

SPECIFICATIONS:

Identity :

Identity

passes test

Ed.: 4 . Vig.: 09.06.2017 .

Prod.: 176131