



PRODUCT CODE: 181040

Nitric Acid 0.1 mol/l (0.1N) volumetric solution

Indicator: Methyl Red

HNO₃

M.= 63,01

CAS [7697-37-2]

EINECS 231-714-2

TARIC 2808 00 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,002 •

BIBLIOGRAPHY:

HAZARDOUS: VLA-EC (HNO₃) 4 ppm10 mg/m³ VLA-ED (HNO₃) 2 ppm5,2 mg/m³

TRANSPORT REGULATIONS: UN: 2031 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l~1,002 kg 1kg~0,998 l

SPECIFICATIONS:

Titer at 20°C
Uncertainty
Traceability

0,999 - 1,001
See certificate
NIST