



**PRODUCT CODE: 181464**

**Silver Nitrate 0.1 mol/l (0.1N) volumetric solution**

Indicator: Potassium Chromate

---

AgNO<sub>3</sub>

M.= 169,87

CAS [7761-88-8]

EINECS 231-853-9

TARIC 2843 21 00 00

**PHYSICAL DATA:** liquid, Clear, Colourless, Miscible with water • D 20/4 1,012 • pH4 - 5 •

**BIBLIOGRAPHY:**

**HAZARDOUS:** C.E: 047-001-00-2 • VLA-ED 0,01 mg/m<sup>3</sup>



H: H412 • H315 • H319 •

P: P273 • P501 • •

**WEIGHT/VOLUME INFORMATION:** 1l~1,012 kg      1kg~0,988 l

**SPECIFICATIONS:**

Titer at 20°C  
Uncertainty  
Traceability

0,999 - 1,001  
See certificate  
NIST

