



PRODUCT CODE: 181535

Potassium Thiocyanate 0.1 mol/l (0.1N) (Reag. USP) volumetric solution

Standardised against Silver Nitrate 0.1M

Indicator: Alum Iron Ammonium

CKNS

KSCN

M.= 97,18

CAS [333-20-0]

EINECS 206-370-1

TARIC 2842 90 80 80

SYNONYMS: Potassium Rhodanide, Potassium Sulphocyanide

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,005 •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral mus 594 mg/kg • LD50 oral rat 854 mg/kg •

WEIGHT/VOLUME INFORMATION: 1l~1,005 kg 1kg~0,995 l

SPECIFICATIONS:

Titer at 20°C
Uncertainty
Traceability

0,999 - 1,001
See certificate
NIST