



PRODUCT CODE: 182112

Nitric Acid 2 mol/l (2N) volumetric solution

Indicator: Methyl Red

HNO₃

HNO₃

M.= 63,01

CAS [7697-37-2]

EINECS 231-714-2

TARIC 2808 00 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,07 • pH1 •

BIBLIOGRAPHY:

HAZARDOUS: C.E: 007-004-00-1 • VLA-EC (HNO₃) 4 ppm10 mg/m³ VLA-ED (HNO₃) 2 ppm5,2 mg/m³



H: H314 •

P: P260 • P264 • P280 • P301+P330+P331 • P303+P361+P353 • P501 • P304+P340 • P305+P351+P338 • P310 • P321 • P338 • P363 • P405 •

TRANSPORT REGULATIONS: UN: 2031 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: F • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l~1,07 kg 1kg~0,93 l

SPECIFICATIONS:

Titer at 20°C
Uncertainty
Traceability

0,999 - 1,001
See certificate
NIST

Ed.: 2 . Vig.: 10.07.2014 .

Prod.: 182112