



PRODUCT CODE: 182115

Silver Nitrate 0.05 mol/l (0.05N) volumetric solution

Indicator: Potassium Chromate

AgNO₃

AgNO₃

M.= 169,87

CAS [7761-88-8]

EINECS 231-853-9

TARIC 2843 21 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,008 •

BIBLIOGRAPHY:

HAZARDOUS: C.E: 047-001-00-2 • VLA-ED 0,01 mg/m³

H: H412 •

P: P273 • P501 •

WEIGHT/VOLUME INFORMATION: 1l~1,008 kg 1kg~0,992 l

SPECIFICATIONS:

Titer at 20°C
Uncertainty
Traceability

0,999 - 1,001
See certificate
NIST