



**PRODUCT CODE: 182116**

**Silver Nitrate 1 mol/l (1N) volumetric solution**

Indicator: Potassium Chromate

AgNO<sub>3</sub>

AgNO<sub>3</sub>

M.= 169,87

CAS [7761-88-8]

EINECS 231-853-9

TARIC 2843 21 00 00

**PHYSICAL DATA:** liquid, Clear, Colourless, Miscible with water • D 20/4 1,138 • pH7 - 9 •

**BIBLIOGRAPHY:**

**HAZARDOUS:** C.E: 047-001-00-2 • VLA-ED 0,01 mg/m3



H: H314 • H411 •

P: P260 • P264 • P273 • P280 • P301+P330+P331 • P501 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P310 • P321 • P338 • P363 • P391 • P405 •

**TRANSPORT REGULATIONS:** UN: 1760 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

**WEIGHT/VOLUME INFORMATION:** 1l-1,138 kg 1kg-0,879 l

**SPECIFICATIONS:**

Titer at 20°C

Uncertainty

Traceability

0,999 - 1,001

See certificate

NIST

