



PRODUCT CODE: 182651

Potassium Permanganate 0.1 mol/l (0.5N) volumetric solution
for the titration of percarbonates

KMnO₄

KMnO₄

M.= 158,04

CAS [7722-64-7]

EINECS 231-760-3

TARIC 2841 61 00 00

SYNONYMS: Permanganic Acid Potassium Salt

PHYSICAL DATA: liquid, Clear, Dark violet, Miscible with water • D 20/4 1,01 •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 1.090 mg/kg •

H: H412 •

P: P273 • P501 •

WEIGHT/VOLUME INFORMATION: 1l~1,01 kg 1kg~0,99 l

OBSERVATIONS: Product controlled as a drug precursor. •

SPECIFICATIONS:

Titer at 20°C

0,999 - 1,001

Uncertainty

See certificate

Traceability

NIST