



PRODUCT CODE: 185560

Silver Nitrate 0.5 mol/l (0.5N) volumetric solution

AgNO₃

AgNO₃

M.= 169,87

CAS [7761-88-8]

EINECS 231-853-9

TARIC 2843 21 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 1,068 •

BIBLIOGRAPHY:

HAZARDOUS: C.E: 047-001-00-2 • VLA-ED (Ag) 0,01 mg/m3



H: H314 • H318 • H411 •

P: P260 • P264 • P273 • P280 • P301+P330+P331 • P501 • P303+P361+P353 • P304+P340 • P305+P351+P338 • P310 • P321 • P338 • P363 • P391 • P405 •

TRANSPORT REGULATIONS: UN: 1760 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l~1,068 kg 1kg~0,936 l

SPECIFICATIONS:

Titer at 20°C

1,000 ±0,001