



PRODUCT CODE: 186983

Silver Nitrate 0.1 mol/l (0.1N) (Reag. USP) volumetric solution

Standardised against Sodium Chloride

Indicator: Eosine

AgNO₃

AgNO₃

M.= 169,87

CAS [7761-88-8]

EINECS 231-853-9

TARIC 2843 21 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,012 • pH4 - 5 •

BIBLIOGRAPHY:

HAZARDOUS: C.E: 047-001-00-2 • VLA-ED 0,01 mg/m3



H: H412 • H315 • H319 •

P: P273 • P280 • P302+P352 • P305+P351+P338 •

WEIGHT/VOLUME INFORMATION: 1l~1,012 kg 1kg~0,988 l

SPECIFICATIONS:

Titer at 20°C

Uncertainty

Traceability

0,999 - 1,001

See certificate

NIST

