


**PRODUCT CODE: 211090**

**2-Propanol 99,7 % for synthesis**

C<sub>3</sub>H<sub>8</sub>O

CH<sub>3</sub>CHOHCH<sub>3</sub>



M.= 60,10

CAS [67-63-0]

EINECS 200-661-7

TARIC 2905 12 00 00

**SYNONYMS:** Dimethylcarbinol, IPA, Isopropanol, iso-Propyl Alcohol, sec-Propyl Alcohol

**PHYSICAL DATA:** liquid, Clear, Colourless, Miscible with water, alcohol, ether and trichloromethane. • Hygroscopic  
 • D 20/20 0,785 • M.P.: -88,5 °C • B.P.: 82,5 °C • n<sub>20</sub>/D : 1,3772 • Flash P.:11,7 °C • Ign. T.:425 °C • Vap. press.  
 (20 °C) 43 hPa • Viscosity 20 °C 2,27 mPa.s • D. M. 20 °C 1,66 Debye • Dielec. constant25 °C 18,3 • Heat evap. 82 °C  
 667 KJ/Kg • Satur. conc. 20 °C 105 g/m<sup>3</sup> • Expl. limit2 %(V)12 %(V) •

**BIBLIOGRAPHY:** Merck Index **13**, 5.228 Sax **INJ000** • Safety **2** , **2952 B** • Kühn-Birett **P 35** • Ullmann **(4.)19** , 433 •  
 Beilstein **1** , **360 I** , **183 II** , **375 III** ,**1439 IV** , **1461** • BRN 635639 • ACS **XI** • ISO 6353/3-1987R - 87 , 62 • BP.**2021** •  
 USP **43** • Ph. Eur. **9.0** (2017) **10.0** (2020) • F.C.C **12 13** • Directive 88/344/E.C.E.92/115/E.C.E.94/52/EC97/60/EC  
 (27/10/1997) 2009/32/CE • Royal Decree 472/1990 (6/4/1990) , 2667/1998 (11 /12/1998), 1101/2011 (22/7/2011), JP  
 XVI •

**HAZARDOUS:** C.E: 603-117-00-0 • RTECS: NT 8050000 • LD50 oral rat >2.000 mg/kg • LC50 rat >20 mg/l / 8h •  
 LD50 skn rbt >2.000 mg/kg • VLA-ED 200 ppm500 mg/m<sup>3</sup> VLA-EC (Spain) 400 ppm1.000 mg/m<sup>3</sup>

**Derived No Effect Level (DNEL)**

Workers Dermal, long exposure (systemic) 888mg/kg Workers Inhalation, long term (systemic) 500 mg/m<sup>3</sup>  
 Population, Dermal, long exposure (systemic) 319mg/kg Population Inhalation, long term (systemic) 89 mg/m<sup>3</sup>  
 Population, oral, long term (systemic) 26mg/kg

**Predicted No Effect Concentration (PNEC)**

Freshwater 140.9mg/l Marine water 140.9mg/l Freshwater Sediment 552mg/kg Marine water Sediment 552mg/kg  
 Soil 28mg/kg



H: H225 • H319 • H336 •

P: P243 • P261 • P271 • P303+P361+P353 • P304+P340 • P305+P351+P338 •

**TRANSPORT REGULATIONS:** UN: 1219 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •  
**WEIGHT/VOLUME INFORMATION:** 1l-0,785 kg 1kg-1,274 l  
**OBSERVATIONS:** Storage below 25°C • Storage away from direct light. •

---

**SPECIFICATIONS:**

Assay (G.C.)	99,7%
Identity :	
Identity	IR passes test
Density at 20/4	0,784-0,786
Acidity	0,0015 meq/g
Alkalinity	0,006 meq/g
Non-volatile matter	0,002 %
Methanol (G.C.)	0,25 %
Water (H2O)	0,2 %