



PRODUCT CODE: 211279

Collodion flexible technical grade

TARIC 3912 20 11 90

Nitrocellulose ~ 4% solution in ethanol-diethyl ether containing ~ 5% castor oil .

PHYSICAL DATA: liquid, Clear, viscous, D 25/4 0,78 • B.P.: 34,6 °C 1.000 hPa • Flash P.: -52 °C • Ign. T.:170 °C • Vap. press. (20 °C) 576 hPa • Expl. limit1,7 %(V)36 %(V) •

BIBLIOGRAPHY: Merck Index **12**, 2.547 13, 2.508 Römp **8** , **2169** • BP.**2020** • USP **42** •

HAZARDOUS:



H: H224 • H336 • H302 • EUH019 • EUH066 •
P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P261 • P264 • P270 • P271 • P280 • P301+P312 • P303+P361+P353 • P304+P340 • P312 • P330 • P370+P378 • P403+P233 • P403+P235 • P405 •

TRANSPORT REGULATIONS: UN: 2059 • ADR: 3/I • IMDG: 3/I • IATA: 3/I • PAX: 351 • CAO: 361 • (B) •

WEIGHT/VOLUME INFORMATION: 1l~0,783 kg 1kg~1,277 l

OBSERVATIONS: Storage between 15 and 25°C •

SPECIFICATIONS:

Density at 25/4

0,775-0,790

Ed.: 5 . Vig.: 07.05.2020 .

Prod.: 211279