



PRODUCT CODE: 211757

Vaseline Soft technical grade
for lubricating glass joints

TARIC 2712 10 90 00

It is a purified mixture of semisolid saturated hydrocarbons having the general formula C_nH_{2n+2} and is obtained from petroleum .

SYNONYMS: White Petrolatum

PHYSICAL DATA: Amorphous mass, White, Insoluble in water M.P.: 38 - 60 °C • n_{60/D} : 1,460 - 1,474 • Flash P.:243 °C •

BIBLIOGRAPHY: BP.2025 • USP 43 •

HAZARDOUS:

SPECIFICATIONS:

Identity :

Identity

Suspension matter

Fatty acids, waxes and resins

Free acids

IR passes test

passes test

passes test

passes test