



PRODUCT CODE: 211841

Magnesium metal, ribbon technical grade

ribbon 3 mm x X0,2 mm

Mg

Mg

M.= 24,31

CAS [7439-95-4]

EINECS 231-104-6

TARIC 8104 90 00 90

PHYSICAL DATA: ribbon, Silvery white, Insoluble in water D 1,738 • M.P.: 651 °C • B.P.: 1.107 °C •

BIBLIOGRAPHY: Merck Index **12**, 5.687 13, 5.675 Sax **MAC750** • Safety **2**, **2177 A** • Römp **8**, **2450** • Kühn-Birett **M 1** • Fieser **1627 2254 3189 4315 5419 6351 7218 10251 11307 12290 13170 16198** •

HAZARDOUS: C.E: 012-002-00-9 • RTECS: OM 2100000 •



H: H250 • H260 •

P: P210 • P223 • P231+P232 • P240 • P241 • P280 • P370+P378 • P402+P404 • P501 •

TRANSPORT REGULATIONS: UN: 1869 • ADR: 4.1/III • IMDG: 4.1/III • IATA: 4.1/III • PAX: 446 • CAO: 448 • (E) •

SPECIFICATIONS:

Assay (Compl.)

99%

Insoluble matter in HCl

passes test

Fe

0,05 %

Ni

0,005 %

Pb

0,005 %

Zn

0,02 %

