



PRODUCT CODE: 216164

Copper(II) Nitrate solution 45% w/w technical grade

CuN_2O_6

$\text{Cu}(\text{NO}_3)_2$

M.= 241,60

CAS [10031-43-3]

EINECS 221-838-5

TARIC 2834 29 40 00

PHYSICAL DATA: liquid, Clear, Blue, Miscible with water • D 20/4 1,513 •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 470 mg/kg •



H: H314 • H410 • H290 •

P: P264 • P270 • P280 • P301+P312 • P302+P352 • P501 • P305+P351+P338 • P321 • P330 • P332+P313 • P337+P313 • P362 •

TRANSPORT REGULATIONS: UN: 3264 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l–1,513 kg 1kg–0,661 l

SPECIFICATIONS:

Assay (as $\text{Cu}(\text{NO}_3)_2$) (Iodom.) 44-46 %

Density at 25/4 $\geq 1,5$