


PRODUCT CODE: 221090

2-Propanol for HPLC gradient / UHPLC supergradient grade

C_3H_8O

$CH_3CHOHCH_3$



M.= 60,10

CAS [67-63-0]

EINECS 200-661-7

TARIC 2905 12 00 00

SYNONYMS: Dimethylcarbinol, IPA, Isopropanol, iso-Propyl Alcohol, sec-Propyl Alcohol

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water, alcohol, ether and trichloromethane. • Hygroscopic
 • D 20/20 0,785 • M.P.: -88,5 °C • B.P.: 82,5 °C • n_{20/D} : 1,3772 • Flash P.:11,7 °C • Ign. T.:425 °C • Vap. press.
 (20 °C) 43 hPa • Viscosity 20 °C 2,27 mPa.s • D. M. 20 °C 1,66 Debye • Dielec. constant 25 °C 18,3 • Heat evap. 82 °C
 667 KJ/Kg • Satur. conc. 20 °C 105 g/m³ • Expl. limit 2 % (V) 12 % (V) •

BIBLIOGRAPHY: Merck Index **13**, 5.228 Sax **INJ000** • Safety **2**, **2952 B** • Kühn-Birett **P 35** • Ullmann **(4.)19**, 433 •
 Beilstein **1**, **360 I**, **183 II**, **375 III**, **1439 IV**, **1461** • BRN 635639 • ACS **XI** • ISO 6353/3-1987R - 87, 62 • BP.**2025** •
 USP **43** • Ph. Eur. **10.0** (2019) **11.0** (2023) • F.C.C **12 13** • Directive
 88/344/E.C.E.92/115/E.C.E.94/52/EC97/60/EC (27/10/1997) 2009/32/CE • Royal Decree 472/1990 (6/4/1990),
 2667/1998 (11/12/1998), 1101/2011 (22/7/2011), JP XVI •

HAZARDOUS: C.E: 603-117-00-0 • RTECS: NT 8050000 • LD₅₀ oral rat >2.000 mg/kg • LC₅₀ rat >20 mg/l / 8h •
 LD₅₀ skn rbt >2.000 mg/kg • VLA-ED 200 ppm500 mg/m³ VLA-EC (Spain) 400 ppm1.000 mg/m³

Derived No Effect Level (DNEL)

Workers Dermal, long exposure (systemic) 888mg/kg Workers Inhalation, long term (systemic) 500 mg/m³
 Population, Dermal, long exposure (systemic) 319mg/kg Population Inhalation, long term (systemic) 89 mg/m³
 Population, oral, long term (systemic) 26mg/kg

Predicted No Effect Concentration (PNEC)

Freshwater 140.9mg/l Marine water 140.9mg/l Freshwater Sediment 552mg/kg Marine water Sediment 552mg/kg
 Soil 28mg/kg



H: H225 • H319 • H336 •

P: P243 • P261 • P271 • P303+P361+P353 • P304+P340 • P305+P351+P338 •

TRANSPORT REGULATIONS: UN: 1219 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •
WEIGHT/VOLUME INFORMATION: 1l-0,785 kg 1kg-1,274 l
OBSERVATIONS: Storage below 25°C • Storage away from direct light. •

SPECIFICATIONS:

Minimum assay (G.C.)	99,9%
Density at 20/4	0,784-0,786
Maximum limit of impurities	
APHA colour	10
Acidity	0,0001 meq/g
Alkalinity	0,0001 meq/g
Non-volatile matter	0,0002 %
Water (H ₂ O)	0,05 %
Suitability for IR spectrometry:	passes test
Gradient at 235 nm	1 mUA
Gradient at 254 nm	1 mUA
UV Spectrum (1cm cell; Ref.: water):	
Transmittance at 207 (Cut off) nm	≥ 10 %
Transmittance at 220 nm	≥ 80 %
Transmittance at 230 nm	≥ 90 %
Transmittance at 250-400 nm	≥ 99 %
Data of interest in HPLC:	
Rohrschneider Polarity	3,9
Eluotropic value E° (Al ₂ O ₃)	0,82
Sol. H ₂ O in solv. at 20°C	miscible

*** Read after purging with nitrogen.**

For critical jobs, purge with nitrogen.

Microfiltered product (0.2 µm) and bottled under nitrogen atmosphere.