



PRODUCT CODE: 251002

Immersion Oil (CE-IVD) for clinical diagnostics

In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

CAS [8001-79-4] EINECS 232-293-8
TARIC 1515 30 90 00

PHYSICAL DATA: liquid, viscous, Colourless to yellowish, Insoluble in water D 20/4 0,96 • B.P.: 313 °C • n_{20/D} : 1,48 • Flash P.:110 °C •

BIBLIOGRAPHY: Safety 2 , 707 B •

HAZARDOUS: RTECS: FJ 1520050 •

WEIGHT/VOLUME INFORMATION: 1l~0,960 kg 1kg~1,042 l

SPECIFICATIONS:

Identity :	
Identity	IR passes test
Density at 20/4	0,955-0,965
Refractive Index n _{20/D}	1,477-1,481

Maximum limit of impurities

Insoluble matter in C ₂ H ₅ OH	passes test
--	-------------