




PRODUCT CODE: 251085

Ethanol 96% v/v (CE-IVD) for clinical diagnostics

In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

C₂H₆O

CH₃CH₂OH



M.= 46,07

CAS [64-17-5]

EINECS 200-578-6

TARIC 2207 10 00 90

SYNONYMS: Ethyl Alcohol

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water and most of the solvents • D 20/4 0,805 • M.P.: -114,1 °C • B.P.: 78,5 °C • n₂₀/D : 1,361 • Flash P.:14 °C • Ign. T.:~425 °C • Vap. press. (20 °C) 57,26 hPa • Viscosity 1,2 mPa.s • Expl. limit 2,5 %(V) 13,5 %(V) •

BIBLIOGRAPHY: Merck Index **12**, 3.806 13, 3.795 Sax **EFU300** • Ullmann **(5.)**9 , 587 • Beilstein **1 , IV , 1289** • BRN 1718733 • ACS **XI** • BP.**2025** • USP **43** • Ph. Eur. **10.0** (2019) **11.0** (2023) • F.C.C **12 13** • Directive 88/344/E.C.E.92/115/E.C.E.94/52/EC97/60/EC (27/10/1997) 2009/32/CE • DAB 10, Royal Decree 472/1990 (6/4/1990) , 2667/1998 (11 /12/1998), 1101/2011 (22/7/2011), F. U. IX , 74, Ph.F. X , 15 •

HAZARDOUS: C.E: 603-002-00-5 • RTECS: KQ 6300000 • LD50 oral rat 6.200 - 15.000 mg/kg / (OECD 401) • LC50 rat >50mg/m³ / (OECD 403) • VLA-ED 1.000 ppm1.920 mg/m³



H: H225 • H319 •

P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P280 • P303+P361+P353 • P370+P378 • P403+P235 • P264 • P305+P351+P338 • P337+P313 •

TRANSPORT REGULATIONS: UN: 1170 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,805 kg 1kg~1,242 l

OBSERVATIONS: May be subject to special tax. • Storage away from direct light, away from sources of ignition and heat. Storage at temperature below 25°C. •

SPECIFICATIONS:

Assay (G.C.) (v/v)	96%
Identity :	
Identity	IR passes test
Density at 20/4	0,804-0,807

Maximum limit of impurities

Non-volatile matter	0,0005 %
Filtered product (1 mm).	