



PRODUCT CODE: 251767

Wright's Eosin-Methylene Blue dye for clinical diagnostics
for hematology, blood smear staining

TARIC 3204 19 00 90

PHYSICAL DATA: powder, Deep-green, Very slightly soluble in water.

BIBLIOGRAPHY: Safety 2 , **3616D** • Conn's IX 105 •

HAZARDOUS:



H: H302 • H319 •

P: P280 • P301+P312 • P330 • P305+P351+P338 • P337+P313 • P501 •

SPECIFICATIONS:

Identity :

Identity

A 1% , 1 cm, λ_{1max} .

A 1% , 1 cm, λ_{2max} .

$\lambda_{of max}$. ABS in CH₃OH

T.L.C.

Maximum limit of impurities

Loss on drying at 110°C

IR passes test

>500

>1000

523 - 525 nm

passes test

15%

Staining suitability

Ed.: 2 . Vig.: 30.05.2014 .

Prod.: 251767