



PRODUCT CODE: 253139

Citrosol (Substitute of Xylene)(CE-IVD) for clinical diagnostics

In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

C₁₀H₁₆

C₁₀H₁₆

M.= 136,24

CAS [5989-27-5]

EINECS 227-813-5

TARIC 2902 19 00 00

PHYSICAL DATA: liquid, Clear, Colourless to yellowish, Immiscible with water. • D 20/4 0,84 • M.P.: -74 °C • B.P.: 176 °C • n_{20/D} : 1,4743 • Flash P.:48 °C • Ign. T.:255 °C • Vap. press. (20 °C) 2,1 hPa • Expl. limit0,7 %(V)6,1 %(V) •

BIBLIOGRAPHY:

HAZARDOUS: C.E: 601-029-00-7 • RTECS: GW 6360000 • LD50 oral rat 4.400 mg/kg •



H: H226 • H315 • H317 • H410 •

P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P261 • P264 • P272 • P273 • P280 • P302+P352 • P303+P361+P353 • P321 • P332+P313 • P333+P313 • P362 • P363 • P370+P378 • P391 • P403+P235 •

TRANSPORT REGULATIONS: UN: 2052 • ADR: 3/III • IMDG: 3/III • IATA: 3/III • PAX: 355 • CAO: 366 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-0,842 kg 1kg-1,188 l

SPECIFICATIONS:

Identity :

Identity

Density at 20/4

Specific rotation [α] 20/D (without dil.)

IR passes test

0,841-0,843

+113 - +122°

Maximum limit of impurities
Water (H₂O)

0,05 %