



PRODUCT CODE: 253296

Copper(II) Sulfate solution d.1.053 for clinical diagnosis
for hematology, blood density

CuO_4S

CuSO_4

M.= 159,60

CAS [7758-98-7]

EINECS 231-847-6

TARIC 2833 25 00 00

PHYSICAL DATA: liquid, Clear, Bluish, Miscible with water D 20/4 1,053 •

BIBLIOGRAPHY: Merck Index **12**, 2.722 13, 2.682

HAZARDOUS: C.E: 029-004-00-0 •



H: H410 • H318 •

P: P273 • P280 • P310 • P305+P351+P338 • P391 • P501 •

TRANSPORT REGULATIONS: UN: 3082 • ADR: 9/III • IMDG: 9/III • IATA: 9/III • PAX: 964 • CAO: 964 • (-) •

WEIGHT/VOLUME INFORMATION: 1l-1,053 kg 1kg-0,949 l

SPECIFICATIONS:

Density at 20/4

1,053±0,001