



PRODUCT CODE: 253681

Eukitt ®, mounting medium for clinical diagnostics

for microscopy, mounting medium (® Registered trade-mark of O. Kindler GmbH)

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, viscous, Colourless, Soluble in esters and ketones. Soluble in aromatic hydrocarbons.
D 20/4 0,95 • M.P.: <0 °C • B.P.:~ 137 °C • n₂₀/D : 1,495 • Flash P.:23 °C(DIN 51794) • Ign. T.:>250 °C • Vap. press. (20 °C) 8 hPa(xylene) • Viscosity 250 - 450 mPa.s(Brookfield) • Expl. limit1 %(V)8 %(V) • D relativa al vapor20 °C3,7(xylene) •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat >2.000 mg/kg • LC50 rat 22 mg/l / 4h • LD L0 skn rbt >2.000 mg/kg •



H: H226 • H332 • H312 • H315 •
P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P261 • P264 • P271 • P280 • P302+P352 • P303+P361+P353 • P304+P340 • P312 • P321 • P322 • P332+P313 • P362 • P363 • P370+P378 • P403+P235 •

TRANSPORT REGULATIONS: UN: 1866 • ADR: 3/III • IMDG: 3/III • IATA: 3/III • PAX: 355 • CAO: 366 • (D/E) •

SPECIFICATIONS:

Identity :
Identity
Refractive Index n₂₀/D

IR passes test
1,493-1,496

