



PRODUCT CODE: 254101

Fixing for fast staining (Panoptic No. 1)(CE-IVD) for clinical diagnostics

for microscopy, blood smear staining or medullary smear staining, In Vitro Diagnostic medical device class A in compliance to the Regulation (EU) 2017/746.

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Pale violet, Miscible with water D 20/4 0,791 • Flash P.:12 °C • Ign. T.:~455 °C • Vap. press. (20 °C) 128 hPa •

BIBLIOGRAPHY:

HAZARDOUS: LD L0 oral man 143 mg/kg • LD50 oral rat 5.628 mg/kg • LC50 rat 64000ppm (V) / 4h • CEE (methanol) 200 ml/m3260 mg/m3 VLA-EC (CH3OH) 250 ppm333 mg/m3 VLA-ED (CH3OH) 200 ppm266 mg/m3



H: H225 • H331 • H311 • H301 • H370 •

P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P260 • P261 • P264 • P270 • P271 • P280 • P301+P310 • P302+P352 • P303+P361+P353 • P304+P340 • P307+P311 • P311 • P312 • P321 • P322 • P330 • P361 • P363 • P370+P378 • P403+P233 • P403+P235 • P405 •

TRANSPORT REGULATIONS: UN: 1992 • ADR: 3(6.1)/II • IMDG: 3(6.1)/II • IATA: 3(6.1)/II • PAX: 352 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,791 kg 1kg~1,264 l

COMPOSITION:

Crystal Violet	2 mg
Methanol s.q.m	1000 ml

SPECIFICATIONS:

Identity :
Identity

passes test

Maximum limit of impurities
Suitability for hematology stain

passes test