



PRODUCT CODE: 255254

DPX, mounting medium fast (toluene base) for clinical diagnostics
for microscopy, mounting medium

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, heavy, Miscible with xylene. Immiscible with ethanol • D 20/4 0,945 • n₂₀/D : 1,52
• Flash P.:12 °C •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 636 mg/kg • LC50 rat 30mg/l / 4h • VLA-ED 50 ml/m3191 mg/m3



H: H225 • H302 • H315 • H361d • H336 • H373 • H304 •
P: P210 • P233 • P240 • P241 • P242 • P501 • P243 • P261 • P264 • P271 • P280 • P302+P352 • P303+P361+P353 •
P304+P340 • P305+P351+P338 • P312 • P321 • P322 • P332+P313 • P337+P313 • P362 • P363 • P370+P378 •
P403+P233 • P403+P235 • P405 •

TRANSPORT REGULATIONS: UN: 1992 • ADR: 3(6.1)/II • IMDG: 3(6.1)/II • IATA: 3(6.1)/II • PAX: 352 • CAO: 364
• (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,945 kg 1kg~1,058 l

SPECIFICATIONS:

Refractive Index n₂₀/D

1,515-1,525

