



PRODUCT CODE: 255598

Histofluid ®, mounting medium for clinical diagnostics

for microscopy, mounting medium (® Registered trade-mark of Paul Marienfeld GmbH&Co.KG)

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, viscous, Clear, Colourless, Insoluble in water Miscible with xylene, acetone, benzene, toluene and chloroform • D (20 °C) 0,95 • B.P.: 137 °C • Flash P.:23 °C • Viscosity 20 °C 250 mPa.s •

BIBLIOGRAPHY:

HAZARDOUS: LD50 oral rat 4.300 mg/kg • LC50 inh rat 5000 ppm / 4h • VLA-EC 661 mg/m3 VLA-EC dermal resorption 150 ppm VLA-ED 100 ppm441 mg/m3



H: H226 • H315 • H319 • H335 • H373 •

P: P210 • P280 • P303+P361+P353 • P305+P351+P338 • P405 • P501 •

TRANSPORT REGULATIONS: UN: 1307 • ADR: 3/III • IMDG: 3/III • IATA: 3/III • PAX: 355 • CAO: 366 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,950 kg 1kg~1,052 l

SPECIFICATIONS:

Refractive Index n20/D

1,493-1,496

