



PRODUCT CODE: 256237

Histofix® Decalcifier 3 for clinical diagnostics

(® Registered trade-mark of Panreac Química S.L.U.) fast decalcifier for fixed tissues,
Contains EDTA and Hydrochloric Acid

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,046 •

BIBLIOGRAPHY:

HAZARDOUS: LC50 rat 3124ppm (V) / 1h • VLA-EC (HCl) 10 ppm15 mg/m3 VLA-ED (HCl) 7,6 mg/m3 VLA-ED (HCl) 5 ppm



H: H290 • H314 • H318 • H335 •

P: P280 • P305+P351+P338 • P303+P361+P353 • P310 • P321 • P405 • P501 •

TRANSPORT REGULATIONS: UN: 1789 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l~1,046 kg 1kg~0,956 l

SPECIFICATIONS:

Suitability: as decalcifier

passes test