



PRODUCT CODE: 256238

Histofix® Decalcifier 2 for clinical diagnostics

(® Registered trade-mark of Panreac Química S.L.U.) medium decalcifier for fixed tissues

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water • D 20/4 1,048 •

BIBLIOGRAPHY:

HAZARDOUS: LC50 rat 3124ppm (V) / 1h • VLA-EC (HCl) 10 ppm15 mg/m3 VLA-ED (HCl) 7,6 mg/m3 VLA-ED (HCl) 5 ppm



H: H290 •

P: P261 • P264 • P271 • P280 • P302+P352 • P501 • P304+P340 • P305+P351+P338 • P312 • P321 • P332+P313 • P337+P313 • P362 • P403+P233 • P405 •

TRANSPORT REGULATIONS: UN: 3264 • ADR: 8/II • IMDG: 8/II • IATA: 8/II • PAX: 851 • CAO: 855 • (E) •

WEIGHT/VOLUME INFORMATION: 1l-1,048 kg 1kg-0,954 l

SPECIFICATIONS:

Suitability: as decalcifier

passes test