

Specification

Buffer Solution pH 9.180

277123

Physical Description:	Liquid
Product Code:	277123
Product Name:	Buffer Solution pH 9.180
Short Description:	usage: These buffer solutions are useful for the calibration of high precision laboratory pH meters
Headline Comment:	* standardized, according to DIN19266
Specifications:	pH (10°C): 9.276 ± 0.01 pH (20°C): 9.225 ± 0.01 pH (25°C): 9.180 ± 0.01 pH (30°C): 9.136 ± 0.01
Hazard pictograms	
WGK:	nwg
Storage:	RT
Signal Word:	Danger
GHS Symbols:	GHS08
H Phrases:	H360FD
P Phrases:	P201 P202 P280 P308+P313 P405 P501

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340

Specification

Buffer Solution pH 9.180

277123

CS:	38221900
	

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode
DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340