



PRODUCT CODE: 281298

Ethanol / Diethyl Ether 1:1 v/v + Phenolphthalein 15 mg/L for acidity determination in olive oil

TARIC 3822 90 00 00

PHYSICAL DATA: liquid, Colourless, Miscible with alcohol • D 20/4 0,76 • Flash P.: <10,0 °C (closed cup) •

BIBLIOGRAPHY: Métodos Analíticos en Alimentaria (Métodos Oficiales de Análisis) Aceites y grasas, Panreac, Métodos Oficiales de Análisis de Aceites y Grasas I, 51 (1993), Ministerio de Agricultura, Pesca y Alimentación, Standard Methods for the Analysis of Oils, Fats and Soaps, II.D.I (1964), IUPAC •

HAZARDOUS: LD L0 oral man 260 mg/kg • LD50 oral rat 1.215 mg/kg • LC50 rat 73000ppm (V) / 2h • VLA-EC 1.540 mg/m3 VLA-EC (diethyl ether) 500 ppm VLA-ED (ethanol) 1.000 ppm 1.910 mg/m3 VLA-ED (diethyl ether) 400 ppm 1.230 mg/m3



H: H225 • H319 • H336 • EUH019 •
P: P210 • P280 • P303+P361+P353 • P305+P351+P338 • P501 • P405 •

TRANSPORT REGULATIONS: UN: 1993 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •
WEIGHT/VOLUME INFORMATION: 1l-0,76 kg 1kg-1,32 l

COMPOSITION:

Phenolphthalein	1,5 mg
Ethyl Ether	50 ml
Ethanol	50 ml

SPECIFICATIONS:

Identity :	
Identity	passes test
Density at 20/4	0,755 - 0,765

