


PRODUCT CODE: 281326

Phenolphthalein solution 0.2% for volumetric analysis

pH indicator 8,0 colourless; 10,0 red violet

$C_{20}H_{14}O_4$

$C_{20}H_{14}O_4$



M.= 318,33

TARIC 3822 90 00 00

SYNONYMS: 3,3-bis-(4-Hydroxyphenyl)-1-(3H)-iso-Benzofuranone, 3,3-bis-(4-Hydroxyphenyl)Phthalide

PHYSICAL DATA: liquid, Clear, Colourless, Miscible with water D 0,903 • Flash P.:24 °C •

BIBLIOGRAPHY:

HAZARDOUS: VLA-ED (ethanol) 1.000 ppm1.910 mg/m3



H: H225 • H319 •

P: P210 • P280 • P305+P351+P338 • P303+P361+P353 • P403+P235 • P501 •

TRANSPORT REGULATIONS: UN: 1170 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 364 • CAO: 353 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l-0,918 kg 1kg-1,089 l

SPECIFICATIONS:

pH transition interval:
colourless
violet red

8,0
10,0