


PRODUCT CODE: 281327

Phenolphthalein solution 1% (Reag. USP, Ph. Eur.) for volumetric analysis
pH indicator 8,0 colourless; 10,0 red violet

C₂₀H₁₄O₄

C₂₀H₁₄O₄



M.= 318,33

TARIC 3822 90 00 00

SYNONYMS: 3,3-bis-(4-Hydroxyphenyl)-1-(3H)-iso-Benzofuranone, 3,3-bis-(4-Hydroxyphenyl)Phthalide

PHYSICAL DATA: liquid, Clear, clear to slightly yellow, Miscible with water D 0,89 • Flash P.:23 °C • Ign. T.:~425 °C •

BIBLIOGRAPHY:

HAZARDOUS: VLA-ED (ethanol) 1.000 ppm1.910 mg/m3



H: H225 • H350 • H341 • H319 •

P: P280 • P303+P361+P353 • P210 • P305+P351+P338 • P405 • P501 •

TRANSPORT REGULATIONS: UN: 1170 • ADR: 3/II • IMDG: 3/II • IATA: 3/II • PAX: 353 • CAO: 364 • (D/E) •

WEIGHT/VOLUME INFORMATION: 1l~0,903 kg 1kg~1,107 l

SPECIFICATIONS:

Identity :	
Identity	passes test
pH transition interval:	
colourless	8,0
violet red	10,0
Sensitivity test	passes test