



PRODUCT CODE: 281366

Alum Iron Ammonium saturated solution for volumetric analysis
indicator for determination of chloride according to Volhard

$H_4FeNO_8S_2 \cdot 12H_2O$
 $NH_4Fe(SO_4)_2 \cdot 12H_2O$

M.= 482,19 CAS [7783-83-7] EINECS 233-382-4
TARIC 2833 30 00 00

SYNONYMS: Iron Alum, Iron(III) Ammonium Sulphate

PHYSICAL DATA: liquid, Yellow-brownish, Miscible with water • D 20/4 1,178 •

BIBLIOGRAPHY:

HAZARDOUS:

WEIGHT/VOLUME INFORMATION: 1l~1,178 kg 1kg~0,849 l

COMPOSITION:

Ammonium Iron(III) Sulphate 12-hydrate	50 g
Sulphuric Acid 96%	0,1 ml
Water (s.q.m)	100 ml